

How to clean a CGM transmitter



First, make sure you have your watertight tester. Properly connecting the tester to the transmitter ensures that fluids do not come in contact with the connector pins inside the transmitter.

Before you prepare to clean your transmitter, ensure you have the tester available.



Gather these supplies

- Mild liquid soap
- Soft-bristled small toothbrush
- Cotton swab



To clean your CGM transmitter

1. Wash your hands thoroughly.
2. Attach the tester to the transmitter by lining up the flat side of the tester with the flat side of the transmitter. Push the tester into the transmitter until the flexible side arms of the tester click into the notches on both sides. When properly connected, the green light flashes 6 times.
3. Hold the tester, saturate a cotton swab in the medical adhesive remover solution and gently rub any adhesive residue until it is fully removed.
4. Rinse the transmitter under room temperature tap water for at least one minute.
5. Prepare a mild liquid soap solution using 1/2 teaspoon of mild liquid soap per litre of room temperature tap water. With the tester still attached, submerge in the mild liquid soap solution and soak for one minute.

For more information, please visit www.medtronic-diabetes.com.au

Always follow the directions for use



Approved sensor insertion locations vary by age. See user guide.

- Container
- Clean lint-free dry cloths

If you will also be removing adhesive residue from your transmitter, you also need a medical adhesive remover such as Zoff®, Unisolve®, Remove®, and cotton swabs. We recommend that you do this after you've removed the transmitter, and before your next sensor insertion.

6. Hold the tester and brush the entire surface using a soft-bristled small toothbrush.
7. Rinse the transmitter under running tap water for at least one minute, and until all visible liquid soap is gone.
8. Dry the transmitter and tester with a clean, dry cloth. Place them on the cloth and air dry completely.
9. Disconnect from the transmitter by gently squeezing the arms of the tester and pulling away from the transmitter.

Information contained herein is not medical advice and should not be used as an alternative to speaking with your doctor. Discuss indications, contraindications, warnings, precautions, potential adverse events and any further information with your healthcare professional.

For more information, please visit www.medtronic-diabetes.com.au

© 2023 Medtronic. All rights reserved. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

2 Alma Road
Macquarie Park NSW 2113 Tel: +61 2 9857 9000 Fax: +61 2 9889 5167
Toll Free: 1800 777 808
medtronic-diabetes.com.au