

InPen™ app Compatibility Guide

Medtronic

InPen™ App version 7.2.0	
Compatible OS versions	
iOS	14, 15, 16, 17, 18
Android	11, 12, 13, 14
Incompatible devices	
Device	Model Number
All tablets	
LG Stylo Series	
Additional notes	
Bluetooth® performance varies due to differences in phone hardware and operating systems. The InPen App is not recommended on Jailbroken or Rooted devices.	

For InPen™ app, all devices are either fully compatible or unverified unless explicitly listed in the table of “Incompatible devices” above.