

# MiniMed™ Mobile App Compatibility Guide

Medtronic

## MiniMed™ Mobile App version 2.8

### Apple

*(unless in the whitelist below, all devices considered "unverified" for iOS 14, 15, 16, 17, 18)*

Device	Model Number	Compatible OS Versions
iPad mini Wi-Fi Only (5th Gen, 2019)	iPad11,1	14, 17
iPad mini Wi-Fi/Cellular (5th Gen, 2019)	iPad11,2	14, 17
iPhone 11	iPhone12,1	14, 15, 16, 17, 18
iPhone 11 Pro	iPhone12,3	14, 15, 16, 17, 18
iPhone 11 Pro Max	iPhone12,5	14, 15, 16, 17, 18
iPhone 12	iPhone13,2	14, 15, 16, 17, 18
iPhone 12 mini	iPhone13,1	14, 15, 17
iPhone 12 Pro	iPhone13,3	14, 15, 16, 17, 18
iPhone 12 Pro Max	iPhone13,4	14, 15, 16, 17, 18
iPhone 13	iPhone14,5	15, 16, 17, 18
iPhone 13 Mini	iPhone14,4	15, 17
iPhone 13 Pro	iPhone14,2	15, 16, 17, 18
iPhone 13 Pro Max	iPhone14,3	15, 16, 17, 18
iPhone 14	iPhone14,7	16, 17, 18
iPhone 14 Plus	iPhone14,8	16, 17, 18
iPhone 14 Pro	iPhone15,2	16, 17, 18
iPhone 14 Pro Max	iPhone15,3	16, 17, 18
iPhone 15	iPhone15,4	17, 18
iPhone 15 Plus	iPhone15,5	17, 18
iPhone 15 Pro	iPhone16,1	17, 18
iPhone 15 Pro Max	iPhone16,2	17, 18
iPhone 6s	iPhone8,1	14, 15
iPhone 6s Plus	iPhone8,2	14, 15

iPhone 7 (China, Japan)	iPhone9,1	14, 15
iPhone 7 (Global)	iPhone9,3	14, 15
iPhone 7 Plus (China, Japan)	iPhone9,2	14, 15
iPhone 7 Plus (Global)	iPhone9,4	14, 15
iPhone 8 (China, Japan)	iPhone10,1	14, 15
iPhone 8 (Global)	iPhone10,4	14, 15
iPhone 8 Plus (China, Japan)	iPhone10,2	14, 15
iPhone 8 Plus (Global)	iPhone10,5	14, 15
iPhone SE 2 (2020)	iPhone12,8	14, 15, 17
iPhone SE 3 (2022)	iPhone14,6	15, 16, 17
iPhone X (China)	iPhone10,3	14, 15
iPhone X (Global)	iPhone10,6	14, 15
iPhone Xr	iPhone11,8	14, 15, 16, 17
iPhone Xs	iPhone11,2	14, 15, 16, 17
iPhone Xs Max	iPhone11,4	14, 15, 16, 17
iPhone Xs Max Global	iPhone11,6	14, 15, 16, 17
iPod touch (7th Gen, 2019)	iPod9,1	14, 15
iPhone 16 Pro	iPhone17,1	18
iPhone 16 Pro Max	iPhone17,2	18
iPhone 16	iPhone17,3	18
iPhone 16 Plus	iPhone17,4	18

Apple Watch		
Apple Watch Series 4	Series 4	7, 8, 9, 10
Apple Watch Series 5	Series 5	7, 8, 9, 10
Apple Watch Series 6	Series 6	7, 8, 9, 10, 11
Apple Watch Series 7	Series 7	8, 9, 10, 11
Apple Watch Series 8	Series 8	9, 10, 11
Apple Watch Series 9	Series 9	10, 11
Apple Watch Series SE (2020)	Series SE (2020)	7, 8, 9, 10
Apple Watch Series SE 2nd Gen (2022)	Series SE 2nd Gen (2022)	7, 8, 9, 10, 11
Apple Watch Ultra	Ultra	9, 10, 11
Apple Watch Ultra 2	Ultra 2	10, 11
Apple Watch Series 10	Series 10	10, 11

## Samsung

*(unless in the whitelist below, all devices considered "unverified" for Android 9, 10, 11, 12, 13, 14)*

Device	Model Number	Compatible OS Versions
Galaxy A10e	SM-S102DL, SM-A102N, SM-A102U, SM-A102U1, SM-A102W	10
Galaxy A20	SCV46, SCV46-j, SCV46-u, SM-S205DL, SM-A205F, SM-A205FN, SM-A205G, SM-A205GN, SM-A205U, SM-A205U1, SM-A205W, SM-A205YN, SM-S205S, SC-02M, SM-A205S	10
Galaxy A20e	SM-A202F, SM-A202K	10
Galaxy A21s	SM-A217F, SM-A217M, SM-A217N, SM-A217F/DS, SM-A217F/DSN, SM-A217M/DS	10
Galaxy A35	SM-A356E, SM-A356E/DS, SM-A356B, SM-A356B/DS, SM-A356U, SM-A356U1	14
Galaxy A50	SM-A505F, SM-A505FM, SM-A505FN, SM-A505G, SM-A505GN, SM-A505GT, SM-A505N, SM-A505U, SM-A505U1, SM-A505W, SM-A505X, SM-A505YN	10, 11
Galaxy A50s	SM-A507F, SM-A5070, SM-A507N, SM-A507FN, SM-S5070	10, 11
Galaxy A51	SM-S515DL, SM-A515F, SM-A515U, SM-A515U1, SM-A515W, SM-A515X, SM-A515F/DS, SM-A515F/DSM, SM-A515F/DSN, SM-A515F/DST, SM-A515F/N	10, 11
Galaxy A51 5G	SM-A5160, SM-A516B, SM-A516N, SM-A516U, SM-A516U1, SM-A516F, SM-A516B/DS, SM-A516F/DSN	10, 11
Galaxy A51 UW	SM-A516V	10, 11
Galaxy A52 5G	SM-A5260, SM-A526B, SM-A526U, SM-A526U1, SM-A526W, SM-A526B/DS	11, 12
Galaxy A52s 5G	SM-A528N	11, 12
Galaxy A54	SM-A546V, SM-A546U, SM-A546U1, SM-A546B, SM-A546B/DS, SM-A546E, SM-A546E/DS, SM-A5460, SM-A546MSM-A546M/DS	14

Galaxy A55	SM-A556V, SM-A556B, SM-A556B/DS, SM-A556E, SM-A556E/DS, SM-A5560	14
Galaxy A7 (2018)	SM-A750F/DS, SM-A750FN/DS, SM-A750G/DS, SM-A750GN/DS	10
Galaxy A71	SM-A715F, SM-A715W, SM-A715X, SM-A715F/DS, SM-A715F/DSM, SM-A715F/DSN	10, 11
Galaxy A71 5G	SM-A7160, SM-A716B, SM-A716F, SM-A716U, SM-A716U1, SM-A716V, SM-A716B/DS, SM-A716F/DS, SM-A716F/DSN	10, 11
Galaxy Note10	SM-N9700, SM-N970F, SM-N970U1, SM-N970U, SM-N970W, SM-N970N, SM-N970X	10, 11, 12
Galaxy Note10 5G	SM-N971N, SM-N971U	10, 11, 12
Galaxy Note10 Lite	SM-N770F, SM-N770F/DS, SM-N770F/DSM	10, 11
Galaxy Note10+	SCV45, SC-01M, SM-N9750, SM-N975C, SM-N975F, SM-N975U, SM-N975U1, SM-N975W, SM-N975D, SM-N975J	10, 11, 12
Galaxy Note10+ 5G	SM-N9760, SM-N976B, SM-N976N, SM-N976Q, SM-N976U, SM-N976V, SM-N976	10, 11, 12
Galaxy Note10+ 5G (China)	SM-N976F	10, 11, 12
Galaxy Note20	SM-N980F, SM-N980F/DS	10, 11, 12, 13
Galaxy Note20 5G	SM-N9810, SM-N981B, SM-N981N, SM-N981W, SM-N981B/DS	10, 11, 12
Galaxy Note20 5G (USA Unlocked)	SM-N981U1	10, 11, 12
Galaxy Note20 5G (USA)	SM-N981U	10, 11, 12

Galaxy Note20 Ultra	SM-N985F, SM-N985F/DS	10, 11, 12, 13
Galaxy Note20 Ultra 5G	SM-N9860, SM-N986B, SM-N986N, SM-N986U, SM-N986U1, SM-N986W, SC-53A, SCG06	10, 11, 12, 13
Galaxy Note9	SC-01L, SCV40, SM-N9600, SM-N960F, SM-N960N, SM-N960U, SM-N960U1, SM-N960W, SM-N960X	10
Galaxy S10	SC-03L, SCV41, SM-G9730, SM-G9738, SM-G973C, SM-G973F, SM-G973N, SM-G973U, SM-G973U1, SM-G973W, SM-G973X	9, 10, 11, 12
Galaxy S10 5G	SM-G977B, SM-G977N, SM-G977P, SM-G977T	9, 10, 11, 12
Galaxy S10 5G (China)	SM-G9770	9, 10, 11, 12
Galaxy S10 5G (USA)	SM-G977U	9, 10, 11, 12
Galaxy S10 Lite	SM-G770F, SM-G770U1, SM-G770F/DS, SM-G770F/DSM	10, 11, 12
Galaxy S10+	SM-G975W, SM-G9750, SM-G9758, SM-G975F, SM-G975N, SM-G975U, SM-G975U1, SC-04L, SCV42, SM-G975X	9, 10, 11, 12
Galaxy S10+ Olympic Games Edition	SC-05L	9, 10, 11, 12
Galaxy S10e	SM-G9700, SM-G9708, SM-G970F, SM-G970N, SM-G970U, SM-G970U1, SM-G970W, SM-G970X	9, 10, 11, 12
Galaxy S20	SM-G980F, SM-G980	10, 11, 12, 13
Galaxy S20 5G	SC-51A, SCG01, SM-G9810, SM-G981B, SM-G981N, SM-G981V, SM-G981W, SM-G981, SM-G981F	10, 11, 12, 13

Galaxy S20 5G (USA Unlocked)	SM-G981U1	10, 11, 12, 13
Galaxy S20 5G (USA)	SM-G981U	10, 11, 12, 13
Galaxy S20 Ultra	SM-G988B	10, 11, 12, 13
Galaxy S20 Ultra 5G	SM-G9880, SM-G988N, SM-G988Q, SM-G988U, SM-G988U1, SM-G988W, SCG03, SM-G988, SM-G988B/DS	10, 11, 12, 13
Galaxy S21 Ultra 5G	SM-G9980, SM-G998B/DS	
Galaxy S20+	SM-G985F	10, 11, 12, 13
Galaxy S20+ 5G	SC-52A, SCG02, SM-G9860, SM-G986B, SM-G986N, SM-G986U, SM-G986U1, SM-G986W, SM-G986, SM-G986F	10, 11, 12, 13
Galaxy S21	SM-G990F, SM-G990F/DS	11, 12, 13, 14
Galaxy S21 5G	SM-G9910, SC-51B, SM-G991B, SM-G991N, SM-G991Q, SM-G991U1, SM-G991W, SM-G991U, SM-G991B/DS	11, 12, 13, 14
Galaxy S21 FE 5G	SM-G9900, SM-G990B, SM-G990E, SM-G990U, SM-G990W, SM-G990B2, SM-G990U1, SM-G990U2, SM-G990W2	11, 12, 13, 14
Galaxy S21+ 5G	SM-G9960, SM-G996B, SM-G996N, SM-G996U1, SM-G996W, SM-G996U, SM-G996B/DS	
Galaxy S22	SM-S901	12, 13, 14

Galaxy S22 5G	SM-S901U1, SM-S9010, SM-S901B, SM-S901E, SM-S901N, SM-S901U, SM-S901W, SM-S901B/DS, SM-S901E/DS	12, 13, 14
Galaxy S22 Ultra	SM-S908	12, 13, 14
Galaxy S22 Ultra 5G	SM-S9080, SM-S908B, SM-S908E, SM-S908N, SM-S908U1, SM-S908W, SM-S908U, SM-S908B/DS, SM-S908E/DS	12, 13, 14
Galaxy S22+	SM-S906	12, 13, 14
Galaxy S22+ 5G	SM-S9060, SM-S906B, SM-S906E, SM-S906N, SM-S906U, SM-S906U1, SM-S906W, SM-S906B/DS, SM-S906E/DS	12, 13, 14
Galaxy S23	SM-S9110, SM-S911B, SM-S911E, SM-S911N, SM-S911U, SM-S911U1, SM-S911W, SM-S911DS, SM-S911B/DS, SM-S911E/DS	13, 14
Galaxy S23 Ultra	SM-S9180, SM-S918B, SM-S918E, SM-S918N, SM-S918U, SM-S918U1, SM-S918W, SM-S918DS, SM-S918B1, SM-S918B/DS, SM-S918E/DS	13, 14
Galaxy S23+	SM-S9160, SM-S916B, SM-S916E, SM-S916N, SM-S916U, SM-S916U1, SM-S916W, SM-S916DS	13, 14
Galaxy S24	SM-S9210, SM-S921B, SM-S921E, SM-S921N, SM-S921U, SM-S921U1, SM-S921W, SM-S921B/DS, SM-S921E/DS	14
Galaxy S24 Ultra	SM-S9280, SM-S928B, SM-S928E, SM-S928N, SM-S928U, SM-S928U1, SM-S928W, SM-S928B/DS, SM-S928E/DS	14
Galaxy S24+	SM-S9260, SM-S926B, SM-S926E, SM-S926N, SM-S926U, SM-S926U1, SM-S926W, SM-S926B/DS, SM-S926E/DS	14
Galaxy S9	SC-02K, SCV38, SM-G9600, SM-G9608, SM-G960F, SM-G960N, SM-G960U, SM-G960U1, SM-G960W, SM-G960, SM-G960X	9, 10

Galaxy S9+	SC-03K, SCV39, SM-G9650, SM-G965F, SM-G965N, SM-G965U, SM-G965U1, SM-G965W, SM-G965X	9, 10
Galaxy Z Flip 4	SM-F7210, SM-F721B, SM-F721U, SM-F721U1, SM-F721W, SM-F721B/DS, SM-F721N	13

## Google

*(unless in the whitelist below, all devices considered "unverified" for Android 9, 10, 11, 12, 13, 14, 15)*

Device	Model Number	Compatible OS Versions
Pixel 3	Pixel 3	9, 10, 11, 12
Pixel 3 XL	Pixel 3 XL	9, 10, 11, 12
Pixel 3a	Pixel 3A	9, 10, 11, 12
Pixel 3a XL	Pixel 3A XL	9, 10, 11, 12
Pixel 4	Pixel 4	10, 11, 12
Pixel 4 XL	Pixel 4 XL	10, 11, 12, 13, 14
Pixel 4a	Pixel 4A	10, 11, 12
Pixel 4a 5G	Pixel 4A (5G)	10, 11, 12, 13, 14
Pixel 5	Pixel 5	11, 12, 13, 14
Pixel 5a	Pixel 5A	11, 12
Pixel 5a 5G	Pixel 5a 5G	11, 12, 13, 14
Pixel 6	G9S9B, GR1YH, Pixel 6, GB7N6	13, 14, 15
Pixel 6 Pro	G8VOU, GLUOG, GF5KQ, Pixel 6 PRO	13, 14, 15
Pixel 6a	GB17L, G9S9B16, GB62Z, GX7AS, Pixel 6A, G1AZG	13, 14
Pixel 7	Pixel 7, GVU6C, GQML3, GO3Z5	15
Pixel 7 Pro	Pixel 7 PRO, GP4BC, GE2AE, GFE4J	15
Pixel 8	PIXEL 8, GKWS6, G9BQD, GA04851-US, GZPFO, GPJ41	15
Pixel 8 Pro	PIXEL 8 PRO, GC3VE, G1MNV	15
Pixel 9	Pixel 9, G2YBB, GUR25, GIB60	15
Pixel 9 Pro	Pixel 9 Pro, GR83Y, GEC77, GWVK6	15
Pixel 9 Pro XL	Pixel 9 Pro XL, GGX8B, GZC4K, GO57S	15



## LG

*(unless in the whitelist below, all devices considered "unverified" for Android 10)*

Device	Model Number	Compatible OS Versions
G7 ThinQ	LG-G710, LM-G710, LM-G710N, LM-G710VM	10
G8 ThinQ	LM-G820, LM-G820N, LM-G820QM, LMG820QM7, LMG820UM0, LMG820UM1, LMG820UM2, LM-G820UMB	10
V40 ThinQ	LM-V405, LM-V409N	10
V50 ThinQ 5G	LM-V450, LM-V500	10
V50 ThinQ 5G (Global)	LM-V500EM	10
V50 ThinQ 5G (Korea)	LM-V500N	10
V50 ThinQ 5G (USA)	LM-V500XM	10
V50S ThinQ	LM-V510N	10
V60	L-51A	10
V60 ThinQ 5G	LM-V600, LM-V600N, A001LG	10

## Huawei

*(unless in the whitelist below, all devices considered "unverified" for Android 9, 10)*

Device	Model Number	Compatible OS Versions
Mate 20	HMA-AL00, HMA-L09, HMA-L29, HMA-LX9, HMA-TL00	9
Mate 20 Pro	LYA-L0C, LYA-L29, LYA-AL00, LYA-AL10, LYA-TL00	9
P30	ELE-AL00, ELE-L04, ELE-L09, ELE-L14, ELE-L29, ELE-L39, ELE-L49, ELE-TL00	9, 10

## Motorola

*(unless in the whitelist below, all devices considered "unverified" for Android 9, 10, 11)*

Device	Model Number	Compatible OS Versions
Moto G30	moto g30, XT2129-1, XT2129-2, PAML0000IN	11

Moto G6	moto g6, XT1925, XT1925-10, XT1925DL	9
Moto G7	XT1962, XT1962-4, moto g7	9
Moto G7 Play	XT1952-4, XT1952-5, moto g7 play	9, 10
Moto G7 Power	XT1955, XT1955-4, moto g7 power	9, 10
Moto G9 Plus	XT2087-1, moto g9 plus	10

## Sony

*(unless in the whitelist below, all devices considered "unverified" for Android 10, 11)*

Device	Model Number	Compatible OS Versions
Xperia 1	802SO, J8110, J8170, J9110, J9180, J9150	10, 11
Xperia 1 II	SOG01, SO-51A	10
Xperia 5	J8210, J8270, J9210, SO-01M, SOV41, J9260	10, 11
Xperia 5 II	A002SO, SO-52A, SOG02, XQ-AS52, XQ-AS62, XQ-AS72	10
Xperia XZ2	H8216, H8266, H8276, H8296	10
Xperia XZ2 Compact	H8314, H8324, SO-05K	10

## Xiaomi

*(unless in the whitelist below, all devices considered "unverified" for Android 9, 10, 11, 12, 13, 14)*

Device	Model Number	Compatible OS Versions
Mi 11 Lite	M2101K9AG, M2101K9AI	13
Mi 6	MCE16, MI 6	9
Poco F4 GT	21121210G	12
Redmi 13C	2311DRN14I, 23100RN82L, 23106RN0DA, 23108RN04Y	13
Redmi 13C 5G	23124RN87G, 23124RN87I, 23124RN87C	13
Redmi Note 10	M2101K7AG, M2101K7AI	11
Redmi Note 9	M2003J15SC, M2003J15SG, M2003J15SS, CX295CIN, CX295VRD	10
Redmi Note 9 Pro	M2003J6B2G, Redmi Note 9 Pro	10
14	23127PN0CC, 23127PN0CG	14
14 Pro	23116PN5BC	14

Note 13 4G	23129RAA4G, 23124RA7EO	14
Note 13 Pro 4G	23117RA68G, 2312FPCA6G	14
Poco C65	2310FPCA4G	13
Poco M6 Pro	23128PC33I	13
Poco X6	23122PCD1G, 23122PCD1I	14
Poco X6 Pro	2311DRK48G, 2311DRK48I	13
Redmi A3	23129RN51X	13
Redmi K70	23113RKC6C	14
Redmi K70 Pro	23117RK66C	14
Redmi K70E	2311DRK48C	13
Redmi Note 13 Pro	2312DRA50C, 2312CRAD3C, 2312DRA50G, 2312DRA50I	14
Redmi Note 13 Pro+	23090RA98C, 23090RA98G, 23090RA98I	13

## OnePlus

*(unless in the whitelist below, all devices considered "unverified" for Android 9, 10, 11, 12, 13.0\*)*

Device	Model Number	Compatible OS Versions
8	IN2010, IN2013, IN2017, IN2019	12
9	LE2110, LE2111, LE2113, LE2115, LE2117	12, 13.0*
11	CPH2447, CPH2449, CPH2451, PHB110	13
10 Pro	NE2210, NE2211, NE2213, NE2215, NE2217	12, 13.0*
10T	CPH2413, CPH2415, CPH2417	12, 13.0*
6T	ONEPLUS A6010, ONEPLUS A6013	9, 10
8 Pro	IN2020, IN2021, IN2023, IN2025	12
8T	KB2000, KB2001, KB2003, KB2005	11, 12
9 Pro	LE2120, LE2121, LE2123, LE2125, LE2127	12, 13.0*

\* Only Android/OxygenOS **13.0** is fully compatible.

Android/OxygenOS **13.1** has been tested but is **incompatible**.

